

12 June 2023

## **Software update for Agilia VP MC, Agilia SP MC & Agilia SP TIVA to v4.1**

To whom it may concern,

Fresenius Kabi is progressively upgrading Agilia VP MC, Agilia SP MC and Agilia SP TIVA pumps to the latest firmware version which is v4.1. This may have already occurred with your Agilia Connect infusion pump due to either routine preventative maintenance or as a repair event on a pump.

The software version is indicated at the top of the screen when you power on the pump.

Software v4.1 brings about some interface enhancements, extra options on the device when programming drug protocols that can be enabled by your pharmacy team such as variable dilution and also allows for compatibility with auto-documentation features when used with an electronic medical record (eMR) system.

The purpose of this letter is to notify you of the upgrade to your pumps.

When a new software version is available for, or applied to your existing Agilia Connect infusion pumps or Vigilant Software Suite system, it is the responsibility of the hospital and its related team to evaluate any workflow or clinical impact as a result of these changes.

This is due to the unique way in which each hospital uses the Agilia Connect infusion pumps and has setup their drug libraries if this feature is used.

The Instructions for Use (IFU's) for the new and older software versions of the Agilia Connect pumps can be found online on our AgiliaASSIST website. Links below:

<https://www.agiliaassist.com.au/vpmc-ifu>

<https://www.agiliaassist.com.au/spmc-ifu>

<https://www.agiliaassist.com.au/sptiva-ifu>

The differences between the old software and the latest software v4.1 can also be found on our AgiliaASSIST website. Link below:

<https://www.agiliaassist.com.au/software>

If you have any questions please contact your local Fresenius Kabi Account Manager.

Yours sincerely,



Carolyn Kirker  
Infusion & Nutrition Systems Product Manager  
Fresenius Kabi Australia Pty Limited & Fresenius Kabi New Zealand Limited  
Tel: +64 (0)21 228 7463  
Carolyn.kirker@fresenius-kabi.com